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www.medevo.de

Bipolar Connection Cable with Pencil

MG004
 Index: F

Last revision:
 2017-01-19

ENGLISH

1 Scope

Bipolar connection cable with pencil Art. No.: 92888-xx
 Maximum rated voltage of accessory:

| Identification | U _{max} |
|--------------------------------------|------------------|
| Bipolar connection cable with pencil | 1.0 kVp |

See also label or catalog.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used.

2 Intended use

Only skilled medical personnel are permitted to use the bipolar connection cables with pencil. The cables must be connected with the appropriate output of the electrosurgical generator and they serve to connect the generator with bipolar electrodes. They are activated by means of a foot-operated switch.

3 Safety notice

WARNING !

See these Instructions for Use, the label or the current product catalog for the maximum rated voltage of the bipolar connection cables with pencil.

If anything is unclear, contact the manufacturer.

Before initial use, new bipolar connection cables with pencil must be cleaned and also sterilized according to a validated procedure (EN ISO 17665) (In this regard see Section 4 "Cleaning and sterilization"). A visual inspection must be performed before each use (see Section 5 "Visual inspection"). It must be ensured that the cable is correctly connected to the generator. In addition, it is necessary to check whether the electrode is inserted firmly in the bipolar connection cable with pencil. It is not permissible to activate the electrosurgical electrode as long as it is in contact with metal objects and/or optics. Throughout the complete procedure, care must be taken that no flammable substances are present in the immediate vicinity, since otherwise a danger of explosion exists.

The high-frequency current used in electrosurgery may interfere with cardiac pacemakers and implanted heart defibrillators, and so affected patients must consult a cardiologist prior to the operation.

4 Cleaning and sterilization

In view of the design, the materials used and the intended use, a maximum limit cannot be defined for the number of cleaning and sterilization cycles that may be possible.

Preparation for cleaning:

Remove the bipolar connection cables from their packaging. Place them in a container/ device provided for cleaning/ sterilization.

Manual Pre-cleaning:

The bipolar connection cables must be disinfected immediately after each use. They should be thoroughly cleaned with a soft brush or synthetic fleece pad then rinsed, since otherwise particles or dried secretions may adhere to them. This could make subsequent cleaning and sterilization difficult or impossible. It must be ensured that areas with difficult access are thoroughly cleaned then rinsed several times. Bipolar connection cables with pencil which

have an integrated suction channel have to be cleaned with the help of a rinsing hose or a rinsing adapter (e.g. Co. Medisave) for at least 60 seconds. Highly alkaline cleaning agents with pH levels above 10 and below 13 have no detrimental influence on the useful life. This manual pre-cleaning has to be done before the automated cleaning.

Automated cleaning:

Only washing machines/ disinfectors with efficiency tested according to EN ISO 15883 may be used. For cleaning, use only commercially available cleaners, which are suitable for cleaning of plastics parts. Thermal disinfection must be performed for at least 5 minutes at a minimum of 90°C [194°F] to a maximum of 95°C [203°F]. The manufacturer's instructions about the dosage of cleaning agent and selection of program must be followed. If chemical disinfection has been applied, a final rinse with distilled or demineralized water must be performed. The bipolar connection cables with pencil must be mounted securely and protected against mechanical damage during machine cleaning/ disinfection. Furthermore the bipolar connection cables with integrated suction channel have to be rinsed with the help of a rinsing hose or a rinsing adapter (e.g. Co. Medisave).

Sterilization:

The bipolar connection cables with pencil are designed exclusively for steam sterilization in autoclaves. The connection cables must be sterilized at a minimum of 134°C [273°F] and a maximum of 138°C [280°F] in saturated steam during a holding time of at least 5 minutes to at most 20 minutes, then dried in vacuum for at least 10 minutes. Sterilization must be performed in accordance with EN ISO 17665 (Sterilization of Medical Devices in Moist Heat). It is not permissible to sterilize the bipolar connection cables with pencil with hot air, EO gas, gamma radiation or plasma.

NOTE: Before use, the bipolar connection cables must be cooled to room temperature.

Limitation of reconditioning:

The useful life of the product will depend on wear, damage and frequency of reconditioning.

5 Visual inspection

Before each use, the insulation of the bipolar connection cables with pencil must be inspected for pressure points or damage. Connection cables with switching function must also be tested for correct functioning. It is not permissible to use bipolar connection cables exhibiting damage or pressure points.

6 Repair and modification

It is not permissible to repair defective bipolar connection cables with pencils. They must be replaced by new bipolar connection cables. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

7 Packaging, storage, transportation, handling

The connection cables must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-

sealed in film. The bipolar connection cables with pencil must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

8 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

9 Disposal

The bipolar connection cables with pencil, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

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11 About this Instruction for Use

Throughout the period of use of the bipolar connection cables with pencil, the Instructions for Use must be kept freely accessible for every user.

A current version may also be found at www.medevo.de

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Reger Medizintechnik GmbH