DEUTSCH

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Arthroscopic Electrodes reusable (up to 5 x)

MG005
Index: VLast Revision:
2018-08-16

Medevo GmbH
Gewerbestraße 10, 78667 Villingendorf
Telefon: 0741/ 175 75 782
Fax: 0741/ 175 75 784
info@medevo.de www.medevo.de

ENGLISH

1 Scope

Monopolar arthroscopic electrodes, article no.: 800-0 to 800-10, 801-0 to 801-10, 802-0 to 802-15, 803-0 to 803-15, 805-0 to 805-15, 806-0 to 806-15, 807-0 to 807-10, 821-0 to 821-10, 832-0 to 832-5, 833-0 to 833-10, 850-0 to 850-6.

Bipolar arthroscopic electrodes, article no.: 870-0 to 870-15, 871-0 to 871-20, 872-0 to 872-10, 873-0 to 873-10, 875-0 to 875-15, 876-0 to 876-10, 879-0 to 879-10, 881-0 to 881-10.

Maximum rated voltage/accessory:

Article No.	U _{max} / P _{max}	Note
800-0 to 800-10, 801-0 to 801-10, 832-1, 832-5, 833-1 to 833-9	1,3 kVp/120 W	Only to be used with REGER cable, article no. 92080, article no. 92083, 92181
802-0 to 802-5, 802-7, 802-10 to 802-13, 803-0 to 803-5, 803-7, 805-0 to 805-5, 805-7, 805-13, 806-0 to 806-5, 806-7, 807-0 to 807-10, 832-0 to 832-5, 833-0 to 833-10, 850-0 to 850-6	1,3 kVp/ 120 W	-
802-6, 802-8, 802-9, 802-15, 803-6, 803-8, 803-9, 805-6, 805-8, 805-9, 805-14, 805-15, 806-6, 806-8, 806-9, 821-0, 821-2	4,3 kVp/ 120 W	-
870-0, 870-2, 870-4, 870-8, 870-12, 871-0, 871-1, 871-4 to 871-20, 872-0 to 872-10, 873-0 to 873-10, 878-0, 879-0 to 879-10, 881-0 to 881-10	1,0 kVp/ 180 W	-
870-1, 870-3, 870-5, 870-6, 870-7 to 870-15, 875-0 to 875-15, 876-0 to 876-10	1,0 kVp/ 120 W	 CAUTION Not to be used for vaporization, only for cutting.

See also label or catalog.

In any combination with another electro-surgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used.

2 Intended use

Only for the use by skilled medical professionals, who have been trained on the usage of the arthroscopic electrodes mentioned in Section 1 "Scope". Read the entire IFU carefully prior to application.

Arthroscopic electrodes are intended for arthroscopic use. They are designed for cutting, vaporizing and coagulating biological tissue (please refer to the table above for exceptions), using 0.9% saline solution. They are not intended for direct contact with the heart or central circulatory system. The bipolar arthroscopic electrodes are indicated for the use in large joints such as knee, shoulder and hip. The arthroscopic electrodes are connected by an electro-surgical cable, electro-surgical handle or handle with suction to the monopolar/ bipolar output of the electro-surgical generator. They are activated by means of a foot-operated switch or electro-surgical handle.

Monopolar arthroscopic electrodes can be operated in combination with the following electro-surgical generators:

ERBE	KLS Martin
Covidien	ValleyLab
EMED	Tekno
BOWA	Berchtold
Medevo ARTro200	

Bipolar arthroscopic electrodes can be operated together with the following electro-surgical generators:

ERBE VIO	KLS Martin Maximum
Mitek VAPR	Medevo ARTro200

The voltage of the used generator may not exceed the maximum rated voltage (See section 1 "Scope") of the arthroscopic electrodes and shall not exceed the frequency 4 MHz. Please follow the instructions for use provided by the manufacturer of the electro-surgical generator.

IMPORTANT: Monopolar arthroscopic electrodes may only be used in combination with a neutral electrode. Prior to usage of the monopolar arthroscopic electrodes it has to be verified that the designated neutral electrode is put on the patient correctly and that the neutral electrode is connected properly with the HF-Generator (electro-surgical generator). The instructions for use of the neutral electrode must be observed.

3 Safety notice - WARNING !

Refer to these Instructions for Use, the label or the current product catalog for the maximum rated voltage of the arthroscopic electrodes.

If anything is unclear, contact the manufacturer.

Prior to each use arthroscopic electrodes must be cleaned, disinfected and also sterilized according to a validated procedure (DIN EN ISO 17665) (In this regard see Section 4 "Cleaning, disinfection and sterilization"). A visual inspection must be performed before each use (see Section 5 "Visual inspection"). It must be ensured that the arthroscopic electrode is inserted firmly in the handle. This must be done carefully, in order to avoid damage to the arthroscopic electrode and/or injuries to the patient or surgical personnel. The arthroscopic electrode may be damaged if excessive force is applied. When using bipolar arthroscopic electrodes, the active electrode and the neutral electrode must be completely surrounded with conductive irrigation fluid (0.9% saline or Ringer's solution).

IMPORTANT: In regards to the application of bipolar arthroscopic electrodes, NEVER use nonconductive irrigation fluids (sterile water, glycine, Purisole, etc.). Always ensure an adequate flow of irrigation fluid (at least 100 mL/min) to the joint.

When applying monopolar arthroscopic electrodes, electrolyte-free solutions, e.g. Purisole, can be used.

It is not permissible to activate the arthroscopic electrode as long as it is in contact with metal objects and/ or optics. During the whole application, make sure the patient will not come into contact with grounded metal objects such as surgical table frames, instrument plates, etc.

Throughout the complete procedure, care must be taken that no flammable substances such as anesthetics, oxidizing gases, endogenous gases etc. are present in the immediate vicinity, since otherwise a danger of explosion exists. Proceed with non-flammable agents for disinfection and cleaning. Do not use alcohol based fluids or similar. Make sure all oxygen connections are dense and leak proof throughout the complete duration of the surgery. After the electro-surgical current has been turned off, the electrode tip may still be hot enough to cause burns.

The high-frequency current used in electro surgery may interfere with cardiac pacemakers and implanted heart defibrillators, and so affected patients must consult a cardiologist prior to the operation.

4 Preparation and application

Make sure the patient is prepared and stored correctly. That means, use insulating table cloth that is dry, absorbent and liquid tight. Insulate the leading areas and contact points towards the patient. It is necessary to use dry cellulose tissues for skin folds, breast folds and between the extremities. Fluid that potentially pooled in body cavities has to be removed prior to the application.

Furthermore it has to be safeguarded that the HF-handle/cable intended for the application has been connected properly to the HF-generator and that the correct performance adjustment has been set. Please follow the IFU of the electro-surgical generator, handle/cable.

Potentially all electro-surgical instruments may cause muscle stimulation during their application. The Medevo arthroscopic electrodes have been designed in such way that the risk of this undesired effect is minimized; however it is still possible that muscle stimulation causes an unexpected movement of the patient in the surgery area.

Do not touch the tip of the electrode during the application. After switching off the HF current the tip of the electrode may still be hot and cause burns. As soon as the arthroscopic electrode is not applied, it shall be deposited on a dry, clean, nonconductive and well visible place that is not in contact with the patient. Any unintended activation of the instrument can lead to burns of the patient.

At completion of the surgery/application the completeness of the system/instrument has to be confirmed.

5 Cleaning, disinfection and sterilization

In view of the design, the materials used and the intended use, the electrodes must not be cleaned, disinfected and sterilized more than up to 5 times. During proper use of the electro-surgical handles, the products also are subject to natural wear and tear depending on type and duration of their application. Therefore a visual inspection (see Section 5 "Visual inspection") has to be done prior to each application.

Medevo recommends automated cleaning and disinfection. A manual cleaning cannot be recommended due to considerably limited efficiency. Do not use sharp pointed objects to clean the product. Take care to thoroughly remove disinfectants and rinse the instruments thoroughly. Medevo recommends the validated treatment processes described in the following. Equivalent deviant processes are possible. It is the sole responsibility of the user to safeguard the suitability of the actually applied procedure by suitable means, e.g. validation, routine examination, verification of material compatibility, etc.

Preparation for cleaning:

Remove the arthroscopic electrodes from their packaging. Place them in a container/ device provided for cleaning/ sterilization.

It is not necessary to dismantle the arthroscopic electrodes.

Manual Pre-cleaning:

Arthroscopic electrodes have to be cleaned immediately after each usage. Clean the arthroscopic electrodes thoroughly with a soft brush or synthetic fleece pad and rinse under running since otherwise particles or dried secretions may adhere to them. This could make subsequent cleaning and sterilization difficult or impossible. It must be ensured that areas with difficult access are thoroughly cleaned then rinsed several times.

Cavities and und lumen (e.g., rinsing/suction channel) have to be rinsed intensively using at least 3 x 20 ml cold mains water (<40°C) with the aids of a rinsing adapter (e.g. from the company Medisafe), or with a syringe or with a water pressure pistol (>30 Sec.).

This pre-cleaning step has to be done prior to the further manual cleaning or prior to cleaning with the cleaning and disinfection machine.

A) Manual cleaning and disinfection:

For cleaning, please prepare an immersion bath with a suitable fluid cleaning agent. Use a disinfection agent that is compatible to the cleaning agent, take care to use a disinfection agent that also can be used for immersion baths. Follow the instructions and recommendations of the manufacturer of the cleaning and disinfection agent. Only use agents that are suitable to clean medical devices made from metal and plastics with a pH-value of between 5.5 and 12.3. It is recommended to use the cleaning agent gigazyme® (Schülke & Mayr).

- According to the instructions of the manufacturer please prepare two separate immersion baths, one for cleaning and one for disinfection using the appropriate and respective agents.
- Put the arthroscopic electrode completely into the cleaner (e.g. 0.5% gigazyme®) in the ultrasonic bath applying an ultrasonic time of 5 Min. and a frequency of 35 KHz. Please follow the instructions of the manufacturer of the cleaning agent. Take care that the instruments do not touch any other parts inside the ultrasonic bath. Take care that sonic shadows are avoided inside the ultrasonic bath.
- Then, clean the arthroscopic electrodes with a soft brush under cold fluent town water (<40°C). Cavities and lumen have to be rinsed intensively by a water pressure pistol (>30 Sec.) or similar. Afterwards rinse the arthroscopic electrodes thoroughly for at least 1 minute under mains water (<40°C) in order to remove any residues of the cleaning agent.
- Inspect the arthroscopic electrodes visually regarding remaining stain. In case soiling/stain is still detectable visually, repeat the aforementioned steps as long as it needs to remove any and all stain.
- Afterwards put the arthroscopic electrodes into the disinfection bath completely. Keep on the contact time set forth by the manufacturer of the disinfection agent, e.g. Korsolex Plus, 3% for 15 min. make sure that all areas and parts of the arthroscopic electrodes are covered by the

disinfectant. Cavities and lumen have to be rinsed with the disinfectant for several times using a syringe e.g. and at least min. 3 times using 20 ml.

- Then, rinse the instruments thoroughly for at least 1 minute with demineralized water in order to remove any and all residues of the disinfectant. Regarding arthroscopic electrodes that have a rinsing/suction channel, please rinse this channel at least 3 times using each 20 ml demineralized water. Rinse the channel with a rinsing adapter (e.g. from the company Medisafe) or with a syringe.
- Dry the arthroscopic electrodes with a lint free wipe and sterile compressed air. Dry the cavities and channels with sterile compressed air as well.

B) Automated cleaning and disinfection

Only apply cleaning and disinfecting machines that are efficiently tested according to DIN EN ISO 15883. Please refer to the information of the manufacturer of the cleaning or disinfecting agent, and only use agents that are suitable to be applied for medical devices made from metal and steel and have a pH-value between 5.5 and 12.3. It is recommended to use neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG). Select the program for thermal disinfection. Follow the instructions for use provided by the manufacturer of the cleaning and disinfecting machine. The arthroscopic electrodes must be mounted securely and protected against mechanical damage during machine cleaning/ disinfection. Do not clean arthroscopic electrodes together with sharp-edged or pointed objects. Put the arthroscopic electrodes in a suitable rinsing device and rinsing adapter if applicable.

Start the program course with following features:

- Pre rinse with cold water for 1 Min.
- Emptying
- Pre rinse with cold water for 3 Min.
- Emptying
- Wash at 55°C for 5 Min. with a 0.5% alkaline cleaner
- Emptying
- Neutralize with warm mains water (>40°C) and Neutralizer (0.1% Neodisher®Z) for 3 Min.
- Emptying
- Intermediate flushing: rinse with warm demineralized water (>40°C) for 2 Min.
- Emptying

Disinfection:

Proceed with the automated thermic disinfection according to the national requirements regarding the AO-value (see ISO 15883, AO-value >3000): 5 Min. for 92°C +/- 2°C.

Drying:

- 30 Min. at 90°C
- Remove rinsing adapter

At end of program course, remove the arthroscopic electrodes and rinsing adapter if applicable, and verify if there are any staining residues. If staining/bonded residues are present, repeat the automated cleaning and disinfection process as long as no visible residues are present. Dry any cavities and insufficiently dried areas with sterile compressed air < 2 bar. After the arthroscopic electrodes are taken out (after additional re-drying on a clean place), immediately pack them in a single-use sterilization package (double package) made of paper/foil or put them in a sterilization container (according to DIN EN ISO 11607 and DIN EN 868).

Sterilization:

Arthroscopic electrodes are designed exclusively for steam sterilization in autoclaves (fractionated fore-vacuum procedure with sufficient product drying). Arthroscopic electrodes have to be sterilized at a minimum of 134°C and maximum of 137°C in saturated steam during a holding time of at least 3 minutes to at most 20 minutes, then dried in a vacuum and sterilized for at least 10 minutes. Regarding the sterilization the Norm DIN EN ISO 17665 (sterilization of medical devices in humid heat) has to be considered. The sterilizer to be used has to fulfill the requirements set forth in DIN EN 13060. A suitable sterilizer is the autoclave 3870 EHS manufactured by Tuttnauer, for example. It is not permissible to sterilize the arthroscopic electrodes with hot air, EO gas, gamma radiation or plasma.

Please follow the recommendations and instructions of the manufacturer regarding loading, handling and drying times of the sterilization device.

NOTE: Before use, the arthroscopic electrodes must be cooled to room temperature.

It is the sole responsibility of the user to safeguard and maintain the sterile condition of the arthroscopic electrodes after the sterilization process.

In case the aforementioned described and recommended chemicals and machines for automated cleaning and disinfection are not available, it is the responsibility of the user to validate their procedure. Also, if another sterilization procedure or sterilizers according to other standards are chosen, the procedure or device deviant from the procedure or device described in section 4 has to be validated and qualified by the user accordingly.

Reconditioning Limitation

The lifetime of the product will depend on natural wear and tear, damage and frequency of reconditioning. The arthroscopic electrodes mentioned in Section 1 "Scope" must not be cleaned, disinfected and sterilized more than up to 5 times maximum. See also paragraph 5 "visual inspection".

6 Visual inspection

Prior to each use, the insulation of the arthroscopic electrodes must be inspected for pressure points or damage. It is not permissible to use arthroscopic electrodes exhibiting damage or pressure points.

At completion of the surgery residues of tissue may adhere to the distal end of the active electrode or sooting may be visible at the distal end. Such adhering tissue or sooting is not regarded as a replacement cause. Additionally, due to the plasma seam and due to mechanical force, deformation of the device may occur or consumption of the insulating material. Also such consumption does not permit a replacement. If it is detected during a long application/surgery that insulating material is consumed, the operator has to exchange the used arthroscopic electrode by a new one. Arthroscopic electrodes that exhibit consumption of insulating material are not allowed to be reprocessed, they have to be exchanged by new ones and have to be disposed of according to paragraph 9 of this IFU „disposal“.

7 Repair and modification

It is not permissible to repair defective arthroscopic electrodes. They must be replaced by new arthroscopic electrodes. Unauthorized modifications and

Arthroscopic Electrodes reusable (up to 5 x)

MG005
Index: V

Last Revision:
2018-08-16

Medevo GmbH
Gewerbstraße 10, 78667 Villingendorf
Telefon: 0741/ 175 75 782
Fax: 0741/ 175 75 784
info@medevo.de www.medevo.de

CE 0483

ENGLISH

repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

A blockage of the suction channel does not permit a replacement.

Arthroscopic electrodes which have a hook as active part are never allowed to be bent.

8 Packaging, storage, transportation, handling

The arthroscopic electrodes must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The arthroscopic electrodes must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

9 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

10 Disposal

The arthroscopic electrodes, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

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12 About this Instructions for Use

Throughout the period of use of the arthroscopic electrodes, the Instructions for Use must be kept freely accessible for every user.

A current version may also be found at www.medevo.de

©2018.Medevo s.r.o.. All rights reserved. All company names, trademarks, service marks, trade names, product names and logos appearing on this document are the property of their respective owners and should be treated accordingly.