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Copyright by Medevo GmbH
All rights reserved
Subject to technical changes without notice

ENGLISH

1 Scope

Bipolar arthroscopic electrode inserts for cutting Art. No.: 890-3S, 890-4S
Bipolar arthroscopic electrode inserts Art. No.: 890-1S, 890-2S, 891-1
Maximum rated voltage of accessory:

Article No.	U _{max} / P _{max}	Comments
890-3S, 890-4S	1.0 kVp/ 50 W	Only suitable for cutting. Only use with Medevo-Exchange-handle 92888.
890-1S, 890-2S, 891-1	1.0 kVp/ 120 W	Only suitable for cutting and coagulating Only use with Medevo - Exchange-handle 92888.

See also label or catalog.
In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used.

2 Authorized use

Only skilled medical personnel are permitted to use the arthroscopic electrode inserts.
The arthroscopic electrode inserts are intended for arthroscopic use. They are designed for cutting (see above table for exceptions), vaporizing and coagulating biological tissue using 0.9% saline solution. They are not intended for direct contact with the heart or central circulatory system. The arthroscopic electrode inserts are connected by an electrosurgical cable, electrosurgical handle or handle with suction to the bipolar output of the electrosurgical generator. They are activated by means of a foot-operated switch or electrosurgical handle. The electrodes are for single use only and must not be resterilized or reused.

3 Safety notice

WARNING !

See these Instructions for Use, the label or the current product catalog for the maximum rated voltage of the arthroscopic electrode inserts.
If anything is unclear, contact the manufacturer.
A visual inspection must be performed before use (see Section 4 "Visual inspection"). It must be ensured that the arthroscopic electrode insert is inserted firmly in the handle. This must be done carefully, in order to avoid damage to the arthroscopic electrode insert and/ or injuries to the patient or surgical personnel. The arthroscopic electrode insert may be damaged if excessive force is applied. During use, the active electrode and the neutral electrode (bipolar application) must be completely surrounded with conductive irrigation fluid (0.9% saline or Ringer's solution). NEVER use nonconductive irrigation fluids (sterile water, glycine, Purisol, etc.). Always ensure an adequate flow of irrigation fluid (at least 100 mL/min) to the joint. It is not permissible to activate the arthroscopic electrode insert as long as it is in contact with metal objects and/ or optics. Throughout the complete procedure, care must be taken that no flammable substances are present in the immediate vicinity, since otherwise a danger of explosion exists. After the electrosurgical current has been turned off, the electrode tip may still be hot enough to cause burns.
The high-frequency current used in electrosurgery may interfere with cardiac pacemakers and implanted heart defibrillators, and so affected patients must consult a cardiologist prior to the operation. Take care to preserve the sterile condition of the electrodes when unpacking them. After unpacking them, the operator must ensure that the sterile condition of the product remains.

4 Visual inspection

Before use, the insulation of the arthroscopic electrode inserts must be inspected for pressure points or damage. It is not permissible to use arthroscopic electrode inserts exhibiting damage or pressure points. In addition, check the packaging for damages and check the expiration date. Arthroscopic electrode inserts with a damaged package or an exceeded expiration date also may not be used.

5 Repair and modification

It is only permissible for the manufacturer to repair defective arthroscopic electrode inserts. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty. Arthroscopic electrode inserts which have a hook as active part, are never allowed to be bent. A blockage of the suction channel does not permit a replacement

6 Packaging, storage, transportation, handling

The arthroscopic electrode inserts must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The arthroscopic electrode inserts must always be handled with the utmost care during transportation and storage. This applies in particular for fine tips and other sensitive areas. The sterility of the product can only be guaranteed as long as the packaging is undamaged. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process. Do not sterilize the arthroscopic electrode inserts.

7 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

8 Disposal

The arthroscopic electrode inserts, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

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10 About this Instruction for Use

Throughout the period of use of the arthroscopic electrode inserts, the Instructions for Use must be kept freely accessible for every user.
A current version may also be found at www.medevo.de