


· Für HF- und RF-Generatoren geeignet
Applicable for HF- and RF-generators

· Autoklavierbar (134° C)
Autoclavable (134° C)



Artikel-Nr. item code	711-3		
Maße dimensions	330 mm		
VPE unit	1	Umax	4,3 KVp
Beschreibung description	Spatelektrode Spatula electrode		



Artikel-Nr. item code	711-5		
Maße dimensions	330 mm		
VPE unit	1	Umax	4,3 KVp
Beschreibung description	L-Hakenelektrode L-hook electrode		



Artikel-Nr. item code	711-4		
Maße dimensions	330 mm		
VPE unit	1	Umax	4,3 KVp
Beschreibung description	J-Hakenelektrode J-hook electrode		



Artikel-Nr. item code	711-7		
Maße dimensions	330 mm		
VPE unit	1	Umax	4,3 KVp
Beschreibung description	Knopfelektrode Button electrode		





· Für HF- und RF-Generatoren geeignet
Applicable for HF- and RF-generators

· Autoklavierbar (134° C)
Autoclavable (134° C)



Ø 4 mm



Artikel-Nr. item code	709-3	710-3
Maße dimensions	330 mm	330 mm
VPE unit	1	
Umax	4,3 kVp	
Beschreibung description	Spatelektrode mit Keramikisolation Spatula electrode with ceramic insulation	




Artikel-Nr. item code	709-5	710-5
Maße dimensions	330 mm	330 mm
VPE unit	1	
Umax	4,3 kVp	
Beschreibung description	L-Hakenelektrode mit Keramikisolation L-hook electrode with ceramic insulation	



Artikel-Nr. item code	709-4	710-4
Maße dimensions	330 mm	330 mm
VPE unit	1	
Umax	4,3 kVp	
Beschreibung description	J-Hakenelektrode mit Keramikisolation J-hook electrode with ceramic insulation	




Artikel-Nr. item code	709-7	710-7
Maße dimensions	330 mm	330 mm
VPE unit	1	
Umax	4,3 kVp	
Beschreibung description	Knopfelektrode mit Keramikisolation Button electrode with ceramic insulation	



330 mm

· Für HF- und RF-Generatoren geeignet
Applicable for HF- and RF-generators

· Autoklavierbar (134° C)
Autoclavable (134° C)



Mit universal Absaug-/ Spülanschluss
With universal suction / irrigation connector



Artikel-Nr. item code	715-3		
Maße dimensions	330 mm		
VPE unit	1	Umax	4,3 kVp
Beschreibung description	Spatelektrode Spatula electrode		



Artikel-Nr. item code	715-4		
Maße dimensions	330 mm		
VPE unit	1	Umax	4,3 kVp
Beschreibung description	L-Hakenelektrode L-hook electrode		



Artikel-Nr. item code	715-5		
Maße dimensions	330 mm		
VPE unit	1	Umax	4,3 kVp
Beschreibung description	J-Hakenelektrode J-hook electrode		

