



Burs and bipolar RF-Electrodes *designed for endoscopy*



EndoBur & EndoVap

Single-Use
and
Reusable

RF-Unit and delicate Handpiece



Diamond Bur designed for controlled abrasion and bone smoothing!

EndoVap

Description	90° Spine Electrode	Front-Tip Spine Electrode	Button Spine Electrode	Hook Spine Electrode
Electrodes				
Tube Diameter [mm]	3,6	3,6	3,6	3,6
Length [mm]	270/320	270/320	270/320	270/320

Electrodes compatible with common RF-Units!


EndoBur

Description	Diamond Round Bur	Diamond Round Bur SideGuard	Round Bur	Round Bur SideGuard
Shaver				
Tube Diameter [mm]	4,0/5,5	4,0/5,5	4,0/5,5	4,0/5,5
Length [mm]	270/320	270/320	270/320	270/320

Burs compatible with common manufacturers' systems!

REGER Shaver System, Handpiece and Foot Switch





EndoBur & EndoVap

REGER Medizintechnik GmbH
Gewerbestraße 10
D-78667 Villingendorf
Tel +49 741 270 698-0
www.reger-med.de
info@reger-med.de



REGER
+ medizintechnik gmbh



Spine EndoBurs & StenoBurs




Abbildung Figure	Kupplung Hub	Kürzel Code	Ø Ø	Länge Length	Pit Bull	Wave Cutter	Synovia Cutter	Barrel Bur SideGuard	Diamond Bur non protected	Diamond Bur SideGuard	RoundBur non protected	Round Bur SideGuard	
	REGER Vimex (Barracuda®)	DRVI	5,50	320 mm					825-204-130r	825-605-174r	825-105-200r	820-105-200r	
			5,50	270 mm					825-604-130r	825-105-174r	825-605-130r	820-605-130r	
			4,00	130 mm	829-202-155r	829-702-200r	829-110-130r	820-604-200r					
			4,00	320 mm						825-104-130r	825-604-200r	825-105-130r	824-605-200r
			4,00	270 mm						825-605-200r	825-104-200r	825-605-200r	829-605-174r
			5,50	130 mm	829-202-155r	829-702-200r	829-110-130r	820-604-200r			829-105-200r		829-105-174r
	Arthrex	DRAR	5,50	320 mm									
			5,50	270 mm									
			4,00	130 mm	823-202-155r	823-702-200r	823-110-130r	820-604-130r					
			4,00	320 mm									
			4,00	270 mm									
			5,50	130 mm	823-202-155r	823-702-200r	823-110-130r	820-604-130r			823-105-200r		823-105-174r
	Smith & Nephew	DRSN	5,50	320 mm									
			5,50	270 mm									
			4,00	130 mm	821-702-200r	821-801-085r	821-110-130r	820-204-130r					
			4,00	320 mm						824-104-174r	824-104-200r	824-604-200r	824-605-174r
			4,00	270 mm									
			5,50	130 mm	821-702-200r	821-801-085r	821-110-130r	820-204-130r			821-105-200r		821-605-200r
	J & J NEXTRA	DRNX	5,50	320 mm									
			5,50	270 mm									
			4,00	130 mm	824-202-155r	824-801-085r	824-110-130r	824-204-130r					
			4,00	320 mm									
			4,00	270 mm									
			5,50	130 mm	824-202-155r	824-801-085r	824-110-130r	824-204-130r			824-105-200r		824-105-174r
	Stryker	DRST	5,50	320 mm					827-605-174r	826-604-200r	826-104-174r	826-605-200r	
			5,50	270 mm					825-104-174r	826-104-200r	826-605-174r	826-105-200r	
			4,00	130 mm	822-202-155r	822-101-200r	822-110-130r	822-204-130r					
			4,00	320 mm						822-604-200r	826-204-130r	820-105-174r	826-605-174r
			4,00	270 mm						822-104-200r	827-605-200r	820-104-174r	826-105-174r
			5,50	130 mm	822-202-155r	822-101-200r	822-110-130r	822-204-130r			822-105-200r		822-605-200r
	Conmed Univalec Vimex (old)	DROR	5,50	320 mm									
			5,50	270 mm									
			4,00	130 mm	827-202-155r	827-801-085r	827-110-130r	820-104-200r					
			4,00	320 mm									
			4,00	270 mm									
			5,50	130 mm	827-202-155r	827-801-085r	827-110-130r	820-104-200r			827-105-200r		827-105-174r

Weitere Varianten auf Anfrage
Additional models upon request
18.10.2023

Varianten:
Options:
Example:

r = Reusable s = Single-Use Sterile

829-102-130r 829-201-130s

Bitte bei Ihrer Bestellung angeben
Please specify with your order

www.reger-med.de

D-78667 Villingendorf | Gewerbestraße 10 | Tel +49 741 270 698-0 | info@reger-med.de

Unsere Shaver Blades und Burs sind für die genannten Systeme passend. Die Bezeichnungen der Kupplungen sind die Handelsmarken der jeweiligen Hersteller. Diese Firmen sind von der REGER GmbH unabhängige Unternehmen. Nicht alle Produkte sind in allen Ländern oder Regionen registriert und zum Verkauf zugelassen. Die Verwendungsangaben können auch je nach Land und Region variieren.
 Produktbilder stellen nur die Referenz. Bitte kontaktieren Sie uns für weitere Informationen. Technische Änderungen und Irrtümer vorbehalten.
 Our Shaver Blades and Burs are suitable for the mentioned systems. The designations of the hubs are the trademarks of the respective manufacturers. The companies are an of REGER GmbH independent company. Not all products are registered and approved for sale in all countries or regions. Indications of use may also vary by country and region.
 Product images are for reference only. Please contact us for further information. Technical changes and errors reserved.



Adapter and handles for EndoVap & StenoVap



Equipment for EndoVap & StenoVap					EndoVap								StenoVap									
Adapter kompatibel zu	Anschluss für	Handgriff mit Schalter	Standard Handstück	Filigranes Handstück	90° Spine Electrode		Front-Tip Spine Electrode		Button Spine Electrode		Hook Spine Electrode		45° Electrode	90° Electrode with Suction	70° Electrode	Hook slim Electrode	Hook Electrode	45° Bent Electrode	45° Electrode	45° Bendable Electrode		
Adapter kompatibel with	connection for	handle with buttons	standard handpiece	delicate handpiece	unsterile, single use														sterile, single use			
					L = 320	L = 270	L = 320	L = 270	L = 320	L = 270	L = 320	L = 270	L = 320	L = 270							L = 150	
	Arthrocare	REGER 6-Pin	92889-6	92888-6	92888-16	889-1	882-1	889-2	882-2	889-6	882-6	889-3	882-3	889-5	882-5	891-11s	891-12s	891-13s	891-14s	891-15s	891-16s	891-17s
	Arthrex		92889-14 (blue plug)																			
	BONSS		92889-14 (blue plug)																			
	Ellman		92889-6																			
	Mitek VUE		92889-6																			
	Mitek VAPR		92889-6																			

Weitere Varianten auf Anfrage
Additional models upon request
25.10.2023

Varianten
Options:
Beispiel:
Example:

Varianten:
Options:
Beispiel:
Example:

Varianten:
Options:
Beispiel:
Example:

Bitte bei Ihrer Bestellung angeben
Please specify with your order

www.reger-med.de

D-78667 Villingendorf | Gewerbestraße 10 | Tel +49 741 270 698-0 | info@reger-med.de

Unsere Elektroden sind für die genannten Systeme passend. Die Bezeichnungen der Anschlüsse sind die Handelsmarken der jeweiligen Hersteller. Diese Firmen sind von der REGER GmbH unabhängige Unternehmen. Nicht alle Produkte sind in allen Ländern oder Regionen registriert und zum Verkauf zugelassen. Die Verwendungsangaben können auch je nach Land und Region variieren.

Produktbilder gelten nur als Referenz. Bitte kontaktieren Sie uns für weitere Informationen. Technische Änderungen und Irrtümer vorbehalten.

Our electrodes are suitable for the mentioned systems. The designations of the connections are the trademarks of the respective manufacturers. The companies are an of REGER GmbH independent company. Not all products are registered and approved for sale in all countries or regions. Indications of use may also vary by country and region. Product images are for reference only. Please contact us for further information. Technical changes and errors reserved.